



Share and access medical imaging data securely.



Application Landscape for Patient Imaging Portal



A Product of



MIRACLE ADVANCE
TECHNOLOGIES CO. LTD.

Integrate. Accelerate. Differentiate.



Version 1.0

Document Control

Release No.	Date	Purpose/Change	Author
V 1.0	Aug 21, 2024	Initial Release	Shahaneaz Shaikat
V 1.1	Sep 8, 2025	Updated	Shahaneaz Shaikat





Table of Contents

1. Business Purpose and Data Confidentiality	3
1.1. Security and Patient Data Management	3
2. On-premise BOM(Bill of Materials)	5
3. User Level Access/Rights with Privileges Mapping	5
4. Architecture and diagram	6
4.1. Mray Upload DICOM and Send Information via API	7
4.2. Patient Portal Site Registration and API Secret Key Generation	8
4.3. Patient Portal Generate New Patients & Orders	8
4.4. Display the Study in the Worklist	9
4.5. API Access for Receiving Mray Link	9
4.6. Logging and Auditing	10
5. Application Landscape	11
5.1. Hospital Workflow	11
5.1.1 Login to the Patient Portal	11
5.1.3 How to configure AD/LDAP for New Site	11
5.1.4 Patient Portal Whiteboard	13
5.1.5 Patient Portal case on the whiteboard	13
5.1.6 whiteboard	13
5.1.6.1 View List	13
5.1.6.2 Filter	14
5.1.6.3 Search	14
5.1.6.4 View Patient Information and Update	15
5.1.6.5 Shareable link or QR code generator buttons on the whiteboard	16
5.1.6.6 Count how many times a link or QR Code has been generated	16
5.2. Settings Options	17
5.2.1 Organization Setup	17
5.2.2 Update organization information	17
• General information	18
• API Settings page	18



• Workflow settings(MRay settings)	19
5.2.3 User List of Patient Portal	19
5.2.4 Adding new Users	20
5.2.5 Update Existing Users	21
5.2.6 Exam Information Page	22
5.2.7 Adding New Exam	22
5.2.8 Update Existing Exam	23
6. Technical Documentation	24
6.1 API Security	24



1. Business Purpose and Data Confidentiality


The primary business purpose of the Patient Portal is to facilitate seamless, secure, and confidential interactions between healthcare providers and patients. It is a platform for managing patient data, generating study orders, and sharing medical images through secure, shareable links or QR codes.

1.1. Security and Patient Data Management

The Patient Portal's critical focus is data confidentiality, especially given the sensitivity of health-related information such as patient records and DICOM files.

The confidentiality measures include:

- ☑ **Sensitive Data Management:** Data, including personally identifiable information (PII) and medical records, is anonymized or de-identified before being loaded into development or non-production environments.



- ☑ **Encryption of data in transit and data at rest:** Data in transit and at rest are properly encrypted and meet the hospital's standards and policies.

HIS_REGISTRATION 1

select hr."REG_ID",hr."FNAME",hr."MNAME",hr."I" Enter a SQL expression to filter results (use Ctrl+Space)

ID	AZ FNAME	AZ MNAME	AZ LNAME	AZ SSN	AZ GENDER
1	q7R2gkdPuNqiz/A+wVEkIA==	[NULL]	p0Ls9QTv3aRwLaOt2WOQwg==	4Dqn1rjYRO4sFyPqf2SzwazhtogUb4+mw02aYa2sr1l=	O8X/HyEaS9vpr
2	ogq8wcdflj10VriJdQRpww==	[NULL]	oEIUMQJz/ugEHy9R5ZUq1g==	PFRNMgX5PIZBa7VY4w3bifU6msXmRL6Tl9GmCHh8WmQ=	ZU9uRz0GL0Ril
3	xcNqKlxCGWRCTD6aVYpoHg==	[NULL]	NnEi92JP3XAuEdiY0sBgWA==	BiqAH092uH5MxCoVFGu4g==	BVSSo45q5cHsi
4	n1rS3eZ3vZrAZhxaLIhqGw==	[NULL]	NDHVwHmsdwPa8ZVak3+yQ==	uIK56hT/uf/84EulhyhBdA==	/52uACZ9/aTTz
5	BfyMhE9LPJT9Jv9MdYGig==	[NULL]	mmcfPjo8Qc9szSlcahhfJw==	WytaNlzsXHlrOpjglwplrw==	zqWy1RQQmdr
6	W50kbq2bfgX7lpVZUoLdeA==	kDkkQFa+yRSPIGptsNHcqA==	dqU1Kg9KAlCg32JPhE1T2w==	sVWPhRudcpX+R0rJhH63Gw==	WBEAh/XhdvK0
7	aw8sEhxO1D6+wTu8MQgO2Q==	[NULL]	J1IagJXNNCIZcRpRZW45Q==	g4jMdygS1nwSt0Fjd85OcvKjDTrPRkXoQvJ4FRg1Wnc=	E41Kfmsjox8tuv
8	cmXntcdMSTrbXa6zNGcGdQ==	[NULL]	hRFkWXQXQb4nwYAmByfIQ==	NYv8SpK2l16CcS/+EFhheea/U2QCIX8lpwMbhiuxWnk=	wnqCAV/aSkLEr
9	16l2LrfsMt+oLlLd1j2WQ==	SMt4FrELH2nbCQrh0t6mXA==	PBdZ8LKIXi5OVS0LNO80/Q==	bxBN+UvtcTngoVRf6o8OiQ==	Stbd1vyTj67lpC
10	x7StylmQ2GvyNQVUUarsg==	iKkb0M0lj0l4mVGrAaLe0A==	PKY5B71fB56EvaYZdCy/cQ==	io8SLcmWlWv2dKjZ+1xGbQ==	UFNuemXi15loz
11			wrOIRvyBEG7+PSdKPLMmag==	[NULL]	Ey1i0lzbLrVltK
12			fB0GRBCx/i99RZszlqykPA==	[NULL]	ti3A37V0dwfzI5
13			Ohprrg6F0XjTKkl910fM8A==	[NULL]	VoVWcM5v6m7
14			4sSCYVrn+BbnsPpl7xYDpA==	[NULL]	UCx6fROjQrOwj
15			NWw8F4hUH+fqPgU2mH6pFg==	[NULL]	A/x3vx6uxV/67f
16	ImRvmxKzDBix24DXEPeA2A==	wAIO+12OSSj0NpUsc2DQA==	+gsme3lUOUjKpNwK7W4hQ==	gulBphGmILhYDSbVcdulJg==	ltN3M9a3b0GBh
17	/Qvcj27dsBoucObuSisdpQ==		pHE3XRdlbXhCd3HV8LemPA==	[NULL]	D920oxOX+T06
18	zmHZxRYRKRQuTaaA660eng==		YXIDTE+GnUk5/S9mqruLMA==	[NULL]	7ILWqyc+ZxtTV
19	0xHT8TAhuLXK/CehW5yRyw==		8fVZ3hq+jGtPhKXKXej9Q==	[NULL]	5cYuo0qAhx8Dh

Refresh Save Cancel Export data 200 200+ 200 row(s) fetched - 0.001s, on 2025-09-05 at 15:43:28

- ☑ **Role-based access control (RBAC):** The Patient Portal adheres to the RBAC security model, which assigns permissions based on user roles within an organization. This approach ensures that only authorized users access specific resources, enhancing security by minimizing unauthorized access and data breaches. RBAC streamlines permission management, facilitates auditing, and helps prevent privilege escalation by aligning access rights with job responsibilities.
- ☑ **API Security**

The Patient Portal API employs multiple security measures to ensure safe and authorized access to its features. API Security measures include:

- **Basic Authentication:** Each API request includes a Basic Auth header, which consists of a **Base64-encoded** string representing the username and password. This method helps prevent unauthorized access to the API and ensures that interactions are secure.

- 
- **API Secret Key:** In addition to Basic Authentication, a Secret Key is required for all API interactions. This key acts as an additional layer of security, ensuring that only requests from authorized sources are accepted. The Secret Key must be included in the header of each API request, providing an additional verification step to secure API access further.

2. On-premise BOM(Bill of Materials)

To successfully deploy and run the Patient Portal, the following hardware components are required:

- **Server-Side Requirements**
 - **On-Premise Server**
 - ☒ Minimum of 32GB RAM
 - ☒ 1TB SSD
 - ☒ quad-core processors.
 - **Dedicated Database Server Running Postgresql(Recommended)**
 - ☒ 32GB RAM
 - ☒ 1TB SSD
 - **Dedicated Backup Storage(Recommended)**
 - ☒ 2TB HDD/SSD(SDD is Recommended)




3. User Level Access/Rights with Privileges Mapping

The Patient Portal implements a Role-Based Access Control (RBAC) system that defines specific access levels based on the user's organizational role. Below is a high-level outline of user access levels and their respective privileges:

User Role	Access Level	Privileges
Super Admin	Full Access	☒ All site access ☒ Site signup approval ☒ Create a new site ☒ View all site information
Site Admin	Single Site Full Access	☒ Manage a single site's settings and configurations ☒ Access to all patient records, studies, and user activity logs ☒ Grant/revoke permissions for other users
API User	API Request	☒ Order Create ☒ Link Request

4. Architecture and diagram

The following application landscape outlines the interaction and data flow between various components within the Patient Portal system, incorporating the exchange of DICOM data, patient information, and API integrations. This landscape enables a seamless workflow from initiating a study in Mray to managing and distributing medical data in the Patient Portal.



4.1. Mray Upload DICOM and Send Information via API

- **Objective:** Initiate the collection of patient data, study information, and DICOM files.
- Upload DICOM files and use the API to transmit site, patient, and study information to the Patient Portal.
- **Actions:**
 - Begin by collecting patient demographics and study details.
 - Prepare DICOM images for upload.
 - Upload DICOM images to Mray's server.
 - Send an API request to the Patient Portal, including:
 - **Site Information:** Clinic or hospital identification
 - **Patient Information:** Name, age, medical history
 - **Study Information:** Modality, study details, physician name

Outcome: DICOM files and related information are securely sent to the Patient Portal.



4.2. Patient Portal Site Registration and API Secret Key Generation

- **Objective:** Register new sites and generate API Secret Keys for secure data exchange.
- **Actions:**
 - Site administrators complete the sign-up process in the Patient Portal.
 - The Patient Portal generates a unique API Secret Key for each new site.
 - The key is provided to the site for future interactions with the Patient Portal.

Outcome: New sites are registered, and secure API keys are generated for data transmission.

4.3. Patient Portal Generate New Patients & Orders

- **Objective:** Automatically generate new patient profiles and study orders based on information received from Mray.
- **Actions:**
 - The Patient Portal creates new patient profiles using the information provided by Mray.
 - A new order is generated, detailing the study and linking it to the patient profile.
 - The system generates a unique Mray Link and QR Code associated with the study.

Outcome: New patients and orders are created and linked to their respective studies. A Mray Link and QR Code are generated.



4.4. Display the Study in the Worklist

- **Objective:** Display the newly created patient profile, study details, and Mray Link/QR Code in the Patient Portal's worklist.
- **Actions:**
 - The system displays the new patient and study order in the worklist.
 - The associated Mray Link and QR Code are available for quick access and sharing.

Outcome: The patient's study is visible in the worklist, with the Mray Link and QR Code available for copying and printing.

4.5. API Access for Receiving Mray Link

- **Objective:** Allow external systems and authorized users to access the Mray Link and study information via the API.
- **Actions:**
 - External systems request the Mray Link via the Patient Portal's API.
 - The API provides access to the Mray Link and QR Code.
 - Every access request is logged, including the number of times the study link has been shared or accessed.

Outcome: Authorized users can securely retrieve and share study links. All interactions are logged for auditing purposes.



4.6. Logging and Auditing

- Objective:
To ensure comprehensive tracking of all system and user activities at the application level for security, compliance, and auditing purposes.
- Actions:
 - Log all significant application-level events, including user creation and management, adjustments to access rights, login and authentication events, and changes to system parameters .
 - Record detailed user activity, including timestamps and specific actions related to accessing and sharing study data.

All logs must be securely retained for a minimum of 90 days.


Outcome: A complete and secure audit trail of all application-level interactions is maintained, ensuring regulatory compliance and enabling effective security monitoring and investigation.

5. Application Landscape

5.1. Hospital Workflow

5.1.1 Login to the Patient Portal


Log in using your login credentials(User Name and Password)




The image shows a login form for the 'Patient Portal'. At the top, there is a logo consisting of a yellow circle with a white plus sign and a black circle, followed by the text 'Patient Portal'. Below the logo, it says 'Login to your Account'. There are two input fields: the first is for the username, labeled with a person icon and the text 'IdapUser'; the second is for the password, labeled with a lock icon and masked with dots. Below these fields is a yellow 'Login' button. At the bottom, there is a small copyright notice: 'Copyright 2025 © Miracle Advance Technologies Co. Ltd.'.

5.1.3 How to configure AD/LDAP for New Site

To configure a user with role-based access control, we have to set up AD/LDAP settings. You have to configure whether this site is on-premises or not. Now, you need to enter the server IP, Port number, and Domain Name for LDAP. Once you've adjusted these settings, click the "**Save**" button to finalize the configuration.




After saving the configuration, you can configure RBAC. To configure, click on the RBAC button. Here, you can create groups and decide who can access the Patient Portal.



5.1.4 Patient Portal Whiteboard

After logging in, you should see the Patient Portal whiteboard.



Patient Name ↑↓	Patient Name (Thai) ↑↓	HN ↑↓	Site ↑↓	Gender ↑↓	DOB ↑↓	SSN ↑↓	Email ↑↓
> UNRUNGSAP SUWIT	สุวิทย์ สุนรุ่งทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
> SAMPHAOLAI PHIRAWAN	พิราวรรณ สำเภาลอย	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
> PRASOETWONG NONGYAO	นงเยาว์ ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
> CHAROENSIN PHANIT	ผาณิต เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
> BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

5.1.5 Patient Portal case on the whiteboard

Clicking on the arrow button lets you review a case in detail.




Patient Name ↑↓	HN ↑↓	Gender ↑↓	DOB ↑↓	SSN ↑↓
▼ Boy Libunao	123456	Male	27-09-1951, 72y	

Accession#	Study	Order Date
20240619_0001	[EXPORT] CHEST PA - 0	19-06-2024 01:16 PM

5.1.6 whiteboard

5.1.6.1 View List

The studies are listed on the whiteboard. Once a study is uploaded through DICOM or API, it will appear there.




The screenshot shows the 'Whiteboard' interface with a table of studies. The table has columns for Patient Name, Patient Name (Thai), HN, Site, Gender, DOB, SSN, and Email. The studies listed are:

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
> UNRUNGSAP SUWIT	สุวิทย์ สุนรุ่งทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
> SAMPHAOLAI PHIRAWAN	พิราวรรณ สำเภาลอย	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
> PRASOETWONG NONGYAO	นงเยาว์ ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
> CHAROENSIN PHANIT	ผาณิต เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
> BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

5.1.6.2 Filter

On the whiteboard, you can filter studies by:


- **Date:** Select a From Date and To Date, or use presets such as Today and Last 7 Days.
- **Site:** Choose a location from the Sites dropdown.
- **Keyword:** Type a patient's name or Hospital Number (HN) into the Search field.



The screenshot shows the 'Whiteboard' interface with the filter section highlighted by a red box. The filter section includes fields for From Date, To Date, Search By, Sites, and a Search field. The table of studies is visible below the filter section.

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
> UNRUNGSAP SUWIT	สุวิทย์ สุนรุ่งทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
> SAMPHAOLAI PHIRAWAN	พิราวรรณ สำเภาลอย	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
> PRASOETWONG NONGYAO	นงเยาว์ ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
> CHAROENSIN PHANIT	ผาณิต เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
> BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

Auto-Refresh: Enable the Auto-Refresh toggle to view new studies as they become available.




The screenshot shows the 'Whiteboard' interface. At the top, there are filters for 'From Date' (03-09-2025), 'To Date' (05-09-2025), 'Search By' (Created On), and 'Sites' (SRIPHAT Imaging Center). Below these are buttons for 'Search', 'Today', 'Yesterday', 'Last 3 Days', and 'Last 7 Days'. On the right, there is an 'Auto Refresh (10s)' toggle switch, which is currently turned off. A red arrow points to this toggle. Below the filters is a table with columns: Patient Name, Patient Name (Thai), HN, Site, Gender, DOB, SSN, and Email. The table contains five rows of patient data. At the bottom, there are pagination controls showing '1' of 10 items.

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
> UNRUNGSAP SUWIT	สุวิทย์ อุ่นรุ่งทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
> SAMPHAOLAI PHIRAWAN	พิรารรณ์ สำเภาลอย	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
> PRASOETWONG NONGYAO	นางสาว ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
> CHAROENSIN PHANIT	ผาณิต เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
> BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

5.1.6.3 Search

The search function enables you to locate patient information using the following criteria: Patient Name, Patient Name (Thai), Site, Gender, Date of Birth (DOB), Social Security Number (SSN), and Email address.



The screenshot shows the 'Whiteboard' interface with the search bar at the top right containing the text '4020175'. A red arrow points to the search bar. The table below shows the search results for the patient with HN 4020175.

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
> UNRUNGSAP SUWIT	สุวิทย์ อุ่นรุ่งทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		

5.1.6.4 View Patient Information and Update

Click on the patient's name to view their details.

Whiteboard

From Date: 03-09-2025 To Date: 05-09-2025 Search By: Created On Sites: SRIPHAT Imaging Center Search: Today Yesterday Last 3 Days Last 7 Days Auto Refresh (10s)

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
UNRUNGSAP SUWIT	สุวิทย์ สุนทรทรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
Accession#	Study	Study Date	Created On	Link Expiry			
20250903MR0108	MRI KNEE JOINT (LEFT) - 0	03-09-2025 12:00 AM	03-09-2025 05:35:33 PM	03-09-2026 (363 days)			
> SAMPHAOLAI PHIRAWAN	พิรารรณ์ สำราญ	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
> PRASOETWONG NONGYAO	นงเยาว์ ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
> CHAROENSIN PHANIT	พาณิชย์ เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
> BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

<< < 1 > >> 10

Once clicked, you can modify the patient's details.

Study Patient Information

HN: 4020175 SSN: SSN Date of Birth: 07-03-1951 Age: 74

First Name: UNRUNGSAP Middle Name: Middle Name Last Name: SUWIT Gender: Male

First Name (Thai): สุวิทย์ Middle Name (Thai): Middle Name Last Name (Thai): สุนทรทรัพย์ Email: Email

Update Close

To save changes, click on the **Update** button.

Study Patient Information

HN: 4020175 SSN: SSN Date of Birth: 07-03-1951 Age: 74

First Name: UNRUNGSAP Middle Name: Middle Name Last Name: SUWIT Gender: Male

First Name (Thai): สุวิทย์ Middle Name (Thai): Middle Name Last Name (Thai): สุนทรทรัพย์ Email: Email

Update Close

5.1.6.5 Shareable link or QR code generator buttons on the whiteboard

On the Patient Portal whiteboard, you have the option to share reports with the patient. You can share reports with patients by:

- Copying a secure link
- Copying, emailing, or printing a QR code
- Retrieving access link via the API

Whiteboard

From Date: 03-09-2025 To Date: 05-09-2025 Search By: Created On Sites: SRIPHAT Imaging Center Q Search Today Yesterday Last 3 Days Last 7 Days Q Search Auto Refresh (10s)

Patient Name	Patient Name (Thai)	HN	Site	Gender	DOB	SSN	Email
UNRUNGSAP SUWIT	สุรินทร์ สุนันททรัพย์	4020175	SRIPHAT Imaging Center (IMC)	Male	07-03-1951, 74y		
Accession#	Study	Study Date	Created On	Link Expiry			
20250903MR0108	MRI KNEE JOINT (LEFT) - 0	03-09-2025 12:00 AM	03-09-2025 05:35:33 PM	03-09-2026 (363 days)			
SAMPHALOI PHIRAWAN	พิราวรรณ สำเภาลอย	2115711	SRIPHAT Imaging Center (IMC)	Female	02-10-1969, 55y		
PRASOETWONG NONGYAO	นงเยาว์ ประเสริฐวงศ์	3142526	SRIPHAT Imaging Center (IMC)	Female	14-04-1957, 68y		
CHAROENSIN PHANIT	ผาณิต เจริญศิลป์	2862153	SRIPHAT Imaging Center (IMC)	Female	18-11-1959, 65y		
BOONKHAM THAMMANOON		3108376	SRIPHAT Imaging Center (IMC)	Male	12-11-1990, 34y		

<< 1 >> 10

5.1.6.6 Count how many times a link or QR Code has been generated

A shareable link and QR code can be generated from here. You can also print the QR code from here. Here, you can see the number of times a link has been clicked.


From Date: 01-06-2024 To Date: 19-08-2024 Q Search Auto Refresh (10s)

Patient Name	HN	Gender	DOB	SSN
Boy Libunao	123456	Male	27-09-1951, 72y	
Accession#	Study	Order Date		
20240619_0001	[EXPORT] CHEST PA - 0	19-06-2024 01:16 PM		

5.2. Settings Options

5.2.1 Organization Setup

This page displays all the existing organizations. You can also update organization information on this page.



The screenshot shows the 'Patient Portal' interface for 'Happy Life Hospital'. The left sidebar contains navigation links: Whiteboard, Patient Info, Settings (highlighted), User List, and Exam Info. The 'Settings' section is expanded, showing 'Organization Setup' as the active option. The main content area is titled 'Organization List' and features a search bar and a 'Clear' button. Below these is a table with the following data:

Org Name	Org Type	Parent OrgName	Hospital Type	Package Name	Org Suffix	Org Prefix	Status	Approved On	Action
Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	PRO			Approved		

The URL at the bottom of the page is <https://mymagingportal.com/#/settings/org-list>.

5.2.2 Update organization information

When the update button is clicked, a pop-up page appears, and you can update the general information, API settings, and workflow here.

The screenshot shows a web form titled "Update Organization Information". At the top, there are three tabs: "General Information", "API Setting", and "Workflow". Three red arrows with numbers 1, 2, and 3 point to these tabs respectively. The form contains several input fields and dropdown menus. On the left side, there are fields for "Org Type*" (set to "Hospital"), "Org Name*" (set to "Happy Life Hospital"), "Tel. No.", "Line ID", "Org Suffix", and "Logo" (with a size of 50*50 and a "+ Choose" button). On the right side, there are fields for "Parent Org*" (set to "Miracle Advance Technology"), "Web URL", "Email" (set to "Hospital2Admin@gmail.com"), and "Org Prefix". Below these is a section titled "Organization Address" with fields for "Country", "Province", "District", "Sub District", "Zip Code", and "Address". At the bottom right, there are "Update" and "Close" buttons.

- General information

The general information page allows you to update the organization's information and address.

This screenshot is similar to the one above, but the "General Information" tab is now selected and highlighted with a yellow underline. The form fields and layout are identical, showing the organization's details and address information. The "Update" and "Close" buttons are visible at the bottom right.

- API Settings page

All incoming and outgoing APIs for a specific organization can be updated here.

Update Organization Information

General Information API Setting Workflow

Submit API Settings

Auth API	Auth API	Base URL	Base URL
Auth Type	Select Auth Type	User Name	UserName
Submit API	Submit API	Password	Password
Submit API Type	Select Submit API Type	Secrete Key	eyJhbGciOiJIUzUxMiIsInR5cCI6IkpXVCJ9.eyJ1aW1laWQlCjB

Receive API Settings

Base URL	Receive Base URL	Receive API	Receive API
User Name	Receive User Name	Receive API Type	Select Receive API Type
Password	Receive Password		

Update Close

- Workflow settings(MRay settings)

All the MRay settings can be added or updated from this page.

Update Organization Information

General Information API Setting Workflow

mRay Setting

mRay Full URL	https://start.mray.app?serverurl={serverurl}&login={login}&studyuid={studyuid}&time={time}&key={key}	Server URL	portal.fav.co.th
Login User	napichead	Pre Shared Key	*****
Link shared to Patient Level	☐		

Update Close

5.2.3 User List of Patient Portal

All the users or roles can be reviewed, added, removed, or updated on this page.

atient Portal

Happy Life Hospital

Settings / User List

User List

+ Add New

Clear

Search Keyword

Full Name	Email	Org Name	Org Type	Parent Org	Hospital Type	Mobile No	Created On	Action
Rad2 User	Rad2CMU@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 23, 2024	
Rad CMU	radcmu@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 23, 2024	
Akash Batash	Akash@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 22, 2024	
Rampaipat Penchome	tang_vet64@yahoo.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	866235584	Aug 21, 2023	
Nougat Radvet	nougatvet@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	848915189	Aug 21, 2023	
Rad 7 LN	pridesys.s.8.1.5.5@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	01719392324	Aug 17, 2023	
Rad 6 L	pridesys.8.155@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	+88001719392324	Aug 17, 2023	
Khandaker Akash	akash.cse.10@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	01911946813	Aug 17, 2023	
Mr. Akash	anfurmiacse18125@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	1670047320	Aug 17, 2023	
Mr. Tester	mrakash@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	98767765432	Jun 22, 2023	

1 2 3 4 5

https://myimagingportal.com/#/settings/user-list

5.2.4 Adding new Users

A new user is added by clicking the "Add New" button.

User List

+ Add New

Clear

Search Keyword

Full Name	Email	Org Name	Org Type	Parent Org	Hospital Type	Mobile No	Created On	Action
Rad2 User	Rad2CMU@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 23, 2024	
Rad CMU	radcmu@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 23, 2024	
Akash Batash	Akash@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging		Feb 22, 2024	
Rampaipat Penchome	tang_vet64@yahoo.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	866235584	Aug 21, 2023	
Nougat Radvet	nougatvet@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	848915189	Aug 21, 2023	
Rad 7 LN	pridesys.s.8.1.5.5@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	01719392324	Aug 17, 2023	
Rad 6 L	pridesys.8.155@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	+88001719392324	Aug 17, 2023	
Khandaker Akash	akash.cse.10@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	01911946813	Aug 17, 2023	
Mr. Akash	anfurmiacse18125@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	1670047320	Aug 17, 2023	
Mr. Tester	mrakash@gmail.com	Happy Life Hospital	Hospital	Miracle Advance Technology	Imaging	98767765432	Jun 22, 2023	

1 2 3 4 5



After clicking, a "Add New User" page appears. Before submitting the form, you must fill in all the information marked with a Red Asterisk(*).

Add User

Org Type*	Hospital	Organization*	Happy Life Hospital
User Role*	Select Role	User Level	Select Level
First Name*	First Name	Last Name*	Last Name
Phone No	Phone No	Email*	Email
Password*	Password	Confirm Password*	Confirm Password
Line Id	Line Id	Medical License	Medical License
Country	Select a country	Province	Select a Province
District	Select a District	Sub-District	Select a Sub-District
Subspecialties	Select Subspecialties	Zip Code	Select a Zip Code

5.2.5 Update Existing Users

Update User


User Information Documents

Org Type*	Hospital	Organization*	Happy Life Hospital	Profile Photo
User Role*	Radiologist	User Level	Select Level	
First Name*	Rad2	Last Name*	User	
Phone No	Phone No.	Email*	Rad2CMU@gmail.com	
Password	*****	Status	Approve	
Medical License	Medical License	Line Id	Line Id	
Country	Select a country	Province	Select a Province	
District	Select a District	Sub-District	Select a Sub-District	
Subspecialties	Select Subspecialties	Zip Code	Select a Zip Code	

Update

5.2.6 Exam Information Page

All the exams can be mapped and configured on the Exam Information page.




The screenshot displays the 'Exam Information' page within the 'Patient Portal' for 'Happy Life Hospital'. The page features a sidebar with navigation options: Whiteboard, Patient Info, Settings, Organization Setup, User List, and Exam Info (highlighted in orange). The main content area shows a table of exams with the following data:

Exam ID	Exam UID	Exam Name	Subspecialty	Action
565		CHEST		
564	UNK	UNKNOWN		
560	XX31	Chest		
581		CT WHOLE ABDOMEN (44503)		
587		chest(41004)(portable)		
586		CHEST(41003)		
588		Chest (Check-Up) 41003		
595		CT BRAIN WITH CONTRAST. (44102)		
594		CT Nasopharynx(44260)		
597		CT UPPER ABDOMEN. (44501)		

At the bottom of the table, there is a pagination control showing '1' of 5 items. The 'Add Exam' button is located in the top right corner of the table area.


5.2.7 Adding New Exam

A new exam can be added. After clicking the "Add Exam" button, you should see a pop-up window where you can enter the name, UID, and Subspecialty for the new exam.



5.2.8 Update Existing Exam

The exam can be updated by clicking on the edit button.





6. Technical Documentation

This section provides detailed API documentation, outlining how to create imaging orders and share links and QR codes through the Patient Portal API. It also covers the security measures implemented to protect these interactions.

Please refer to the document for API documentation on creating imaging orders and sharing links and QR codes – [!\[\]\(4729e517bc6a7cd81c8025b9646574fb_img.jpg\) **Patient Portal API Documentation**](#)

6.1 API Security

The Patient Portal API is secured using Basic Authentication (Basic Auth), which requires API users to provide a valid username and password combination to access API functions. This security measure ensures that only authorized users can execute API calls, providing a simple yet effective layer of security.

- **Basic Authentication:** Each API request includes a Basic Auth header. This header consists of a **Base64**-encoded string representing the username and password. This method helps prevent unauthorized access to the API and ensures that interactions are secure.